

# How to Determine System Dwell Volume: Theory and Practice



## ABSTRACT

Determining the dwell volume of an LC system is a useful activity. This technical note describes the theory behind dwell volume and how to measure it accurately for use in gradient method transfers between different instruments to maintain peak retention, selectivity and resolution.

## INTRODUCTION

The dwell volume or gradient delay volume is described as the system volume between the point at which eluents are combined in the pump mixer and the head of the column.

Dwell volumes differ considerably depending on the configuration of the instrument, the vendor, the volume of tubing, between binary or quaternary pumps and if any mixers are installed. Typical instrument dwell values can be seen in Table 1.

Instrument	Typical Dwell Volume (µL)
Agilent HP1100 Binary	180-900
Agilent HP1100 Quaternary	800-1100
Agilent 1200 RRLC	~300
Dionex P680A Quaternary	<400
Thermoquest P4000 Quat	<600
Waters Alliance 2695 Quat	600
Waters Varian 9012 Ternary	1000
Waters Acquity UPLC	~100

**Table 1** Typical dwell volume values for different instrument manufacturers

Dwell volume considerations are less important for isocratic analyses as the mobile phase composition remains constant.

For gradient chromatography it is different. The volume between initial eluent mixing and the point at which the mixed eluent reaches the column can have significant impact on selectivity and retention. Moving gradient methods between instruments can easily result in different chromatography (retention times, peak elution order) as a result of dwell volume differences.

Understanding the dwell volumes for your instruments is therefore good practice to ensure seamless gradient method transfers between instruments.

## METHODOLOGY


The dwell volume may be determined as follows:

1. Replace the column with a Zero Dead Volume (ZDV) connector.
2. Place water on Solvent Line A and water containing 0.1% v/v acetone on Solvent Line B.
3. Monitor the run using UV detection at 265 nm.
4. Run the gradient program in Table 2 at 2 mL/min.

Time	%B
0	0
10	100
12	100
12.5	0
Post: 3 mins	

**Table 2** Dwell volume gradient

The resulting chromatogram is seen in Figure 1. To identify the midpoint of the gradient ( $t_{0.5}$ ), subtract the UV absorbance at the end of the gradient ( $A_{max}$ ) from the UV absorbance at the beginning of the



**Figure 1** Example of chromatogram produced to calculate dwell volume

## CONCLUSION

A simple, generally applicable method for determining dwell volume in gradient chromatography has been described.

gradient ( $A_{min}$ ). This is divided by two and the time of  $A_{0.5}$  is determined. The dwell time is calculated as:

$$t_D = t_{0.5} - \left(\frac{t_G}{2}\right)$$

where  $t_D$  is dwell time (min) and  $t_G$  is gradient time (min). Dwell volume is then calculated as:

$$V_D = t_D \times F$$

where F is the flow rate (mL/min)

System Dwell Volume	
Gradient endpoint ( $A_{max}$ )	123.138 mAU
Gradient midpoint ( $A_{0.5}$ )	61.569 mAU
Gradient midpoint ( $t_{0.5}$ )	5.566 min
Dwell time $t_D = t_{0.5} - \left(\frac{t_G}{2}\right)$	0.566 min
Dwell volume $V_D = t_D \times F$	1.13 mL

**Table 3** Worked example to calculate dwell volume