

CUSTOM PACKAGING & KITTING SOLUTIONS



Custom Specimen Collection Kits

Customer Need: A clinical laboratory needed a kit designed to retrieve specific biological specimens from test subjects in the field and ship directly to labs for analysis.

VWRCATALYST Custom Packaging & Kitting Solutions delivered a kit that included the following materials:


- 3 10.0mL K2 EDTA Vacutainer Tube
- 4 10.0mL Sodium Heparin Vacutainer Tube
- 2 2.5mL PAXgene RNA Vacutainer Tube
- 1 5.0mL SST Vacutainer Tube
- 1 15mL Transfer Tube, Conical, PP, graduated, sterile
- 1 21g Butterfly Safety-Lok Needle
- 1 Standard Needle / Vacutainer Holder
- 1 2" x 2" Non-Sterile Gauze
- 1 Alcohol Prep Pad
- 1 Tourniquet, 1" x 18", Latex Free (Blue)
- 11 Cryolabels, Barcode BLM
- 1 Bag
- 1 95kPa Biohazard Bag
- 1 4-Bay Absorbent AQUI-Pak
- 1 6-Bay Absorbent AQUI-Pak
- 1 Requisition Form, 1 page
- 1 Box for Mini EPS
- 1 UN3373 Category B Label
- 1 Protocol Label
- 1 FedEx Return Service Label

VALUE DELIVERED

- Custom branding
- Full cGMP traceability of lot controlled components
- Packaging scaled to minimize freight costs

UNIQUE FEATURE

- Kit-Track ID label offers comprehensive disposition & tracking information
- On-demand printing of custom requisition forms
- Online order portal for clinical and commercial kit distribution management



For more information or to discuss your unique needs, please e-mail kittingservices.eu@vwr.com

